



# Recombinant Human Nuclear protein localization protein 4 homolog (NPLOC4)

<b>Product Code</b>	CSB-BP844706HU
<b>Abbreviation</b>	NPLOC4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8TAT6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AESIIIRVQ SPDGVKRLTA TKRETAATFL KKVAKKEFGFQ NNGFSVYINR NKTGEITASS NKSLNLLKIK HGDLLFLFPS SLAGPSSEME TSVPPGFKVF GAPNVVEDEI DQYLSKQDGK IYRSRDPQLC RHGPLGKCVH CVPLEPFDED YLNHLEPPVK HMSFHAYIRK LTGGADKGF VALENISCKI KSGCEGHLPW PNGICTKCQP SAILNRQKY RHVDNIMFEN HTVADRFLDF WRKTGNQHFG YLYGRYTEHK DIPLGIRAEV AAIYEPPIQIG TQNSLELLED PKAEVVDEIA AKLGLRKVGW IFTDLVSEDV RKGTVRYSRN KDTYFLSSEE CITAGDFQNK HPNMCRLSPD GHFGSKFVTA VATGGPDNQV HFEGYQVSNQ CMALVRDECL LPCKDAPELG YAKESSESEQY VPDVFKDQVD KFGNEITQLA RPLPVEYLII DITTFPKDP VYTFSSISQNP FPIENRDVLG ETQDFHSLAT YLSQNTSSVF LDTISDFHLL LFLVTNEVMP LQDSISLLE AVRTRNEELA QTWKRSEQWA TIEQLCSTVG GQLPGLHEYG AVGGSTHTAT AAMWACQHCT FMNQPVTGHC EMCSLPRT
<b>Source</b>	Baculovirus
<b>Target Names</b>	NPLOC4
<b>Protein Names</b>	Recommended name: Nuclear protein localization protein 4 homolog Short name= Protein NPL4
<b>Expression Region</b>	2-608
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.