



Recombinant Human F-box only protein 30 (FBXO30)

Product Code	CSB-MP844709HU
Abbreviation	FBXO30
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TB52
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEEELQHS HC VNCVSRRCMT RPEPGISCDL IGCPLVCGAV FHSCKADEHR LLCPFERVPC LNSDFGCPFT MARNKVAEHL EMCPASVVCC TMEWNRWPVS YADRKSYENL SRDVDEVAQL DMALALQDQR MLLESLKVAT MMSKATDKVS KPREQISVKS SVPEIPHANG LSVVDEESYG ALYQATVETT RSLAAALDIL NTATRDIGML NTSVPNDMDE QQNARESLED QNLKDQDHLY EEEIGAVGGI DYNDTNQNAQ SEQNGSSDLL CDLNTSSYDT SALCNGFPLE NICTQVIDQN QNLHGDSKQS NLTNGDCVAS SDGTSKPSSS LAVAAQLREI IPSSALPNGT VQHILMPDDE GEGELCWKKV DLGDVKNVDV LFSHAPSFN FLSNSCWSKP KEDKAVDTS D LEVAEDPMGL QGIDLITAAL LFCLGDSPGG RGISDSRMAD IYHIDVGTQT FSLPSAILAT STMVGEIASA SACDHANPQL SNPSPFQTLG LDLVLECVAR YQPKQRSMFT FVCGQLFRRK EFSSHFKNVH GDIHAGLNGW MEQRCPLAYY GCTYSQRRFC PSIQGAKIIH DRHLRSFGVQ PCVSTVLVEP ARNCVLGLHN DHLSSLPFEV LQHIAGFLDG FSLCQLSCVS KLMRDVCGSL LQSRGMVILQ WGKRKYPEGN SSWQIKEKVV RFSTAFCSVN EWKFADILSM ADHLKKCSYN VVEKREEAIP LPCMCVTREL TKEGRSLRSV LKPVL
Source	Mammalian cell
Target Names	FBXO30
Protein Names	Recommended name: F-box only protein 30
Expression Region	1-745
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called



SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. This protein belongs to the Fbxs class and it is upregulated in nasopharyngeal carcinoma.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.