



Recombinant Human Cytoplasmic dynein 2 light intermediate chain 1 (DYNC2LI1)

Product Code	CSB-MP844722HU
Abbreviation	DYNC2LI1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TCX1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPSETLWEIAKAEVEKRGINGSEGDAEIAEKVFFFIGSKNGGKTTIILRCLDRD EPPKP TLALEYTYGRRRAKGHNTPKDIAHFWELGGGTSLLDLISIPITGDTLRTFSLVVL DLSKP NDLWPTMENLLQATKSHVDKVIKLGKTNAAVSEMRQKIWNNMPKDHPQD HELIDPPV PLVIIGSKYDVFQDFESEKRKVICKTLRFVAHYYGASLMFTSKSEALLLKIRGVIN QLAF GIDKSKSICVDQNKPLFITAGLDSFGQIGSPPVPENDIGKLHAHSPMELWKKVY EKLFP KSINTLKDIDKPARDPQYAENEVDEMRIQKDLELEQYKRSSSSKSWKQIELD
Source	Mammalian cell
Target Names	DYNC2LI1
Protein Names	Recommended name: Cytoplasmic dynein 2 light intermediate chain 1 Short name= Dynein 2 light intermediate chain
Expression Region	1-352aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein of Isoform 2
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.