



# Recombinant Human Putative E3 ubiquitin-protein ligase SH3RF2 (SH3RF2)

<b>Product Code</b>	CSB-BP844731HU
<b>Abbreviation</b>	SH3RF2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8TEC5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDDLTLDDL ECPVCFEKLD VTAKVLPCQH TFCKPCLQRV FKAHKELRCP ECRTPVFSNI EALPANLLLV RLLDGVRSGQ SSGRGGSFRR PGTMTLQDGR KSRTNPRRLQ ASPFRLVPNV RIHMDGVPR KALCNYRGQN PGDLRFNKG IILLRRQLDE NWFYQGEINGI SGNFPASSVE VIKQLPQPPP LCRALYNFDL RGKDKSENQD CLTFLKDDII TVISRVDENW AEGKLGDKVG IFPILFVEPN LTARHILLEKN KGRQSSRTKN LSLVSSSSRG NTSTLRRGPG SRRKVPQGFS ITTALNTLNR MVHSPSGRHM VEISTPVLIS SSNPSVITQP MEKADVPSSC VGQVSTYHPA PVSPGHSTAV VSLPGSQQHL SANMFVALHS YSAHGPDELD LQKGEGVRVL GKCQDGWLRG VSLVTGRVGI FPNNYVIPIF RKTSSFPDSR SPGLYTTWTL STSSVSSQGS ISEGDPQRSR PFKSVFVPTA IVNPVRSTAG PGTLGQGLR KGRSSMRKNG SLQRPLQSGI PTLVVGSLRR SPTMVLRPQQ FQFYQPQGIP SSPSAVVVEM GSKPALTGEP ALTCISRGSE AWIHSAASSL IMEDKEIPIK SEPLPKPPAS APPSILVKPE NSRNGIEKQV KTVRFQNYSP PPTKHYTSH TSGKPEQPAT LKASQPEAAS LGPEMTVLFA HRSGCHSGQQ TDLRRKSALG KATTLVSTAS GTQTVFPSK
<b>Source</b>	Baculovirus
<b>Target Names</b>	SH3RF2
<b>Protein Names</b>	Recommended name: Putative E3 ubiquitin-protein ligase SH3RF2 EC= 6.3.2.- Alternative name(s): Heart protein phosphatase 1-binding protein Short name= HEPP1 Protein phosphatase 1 regulatory subunit 39 RING finger protein 158
<b>Expression Region</b>	1-729
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.