



Recombinant Human PDZ domain-containing protein GIPC3 (GIPC3)

Product Code	CSB-MP844739HU
Abbreviation	GIPC3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TF64
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEGAAAREAR GTETPRASAP PPAPSEPPAA PRARPRLVFR TQLAHGSPTG KIEGFTNVRE LYAKIAEAFG IAPTEILFCT LNSHKVDMQK LLGGQIGLED FIFAHVRGET KEVEVTKTED ALGLTITDNG AGYAFIKRIK EGSINRIEA VCVGDSIEAI NDHSIVGCRH YEVAKMLREL PKSQPFTLRL VQPKRAFDMI GQRSRSSKCP VEAQVTSGRE TLRLRSGGAA TVEEAPSEFE EEASRKVDDL LESYMGIRDP ELASTMVETS KKTASAQEFA RCLDSVLGEF AFPDEFVVEV WAAIGEAREA CG
Source	Mammalian cell
Target Names	GIPC3
Protein Names	Recommended name: PDZ domain-containing protein GIPC3
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.