



Recombinant Bovine coronavirus Protein I (N)

Product Code	CSB-BP844974BJF
Abbreviation	N
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8V431
Product Type	Recombinant Protein
Immunogen Species	Bovine coronavirus (strain 98TXSF-110-LUN) (BCoV-LUN) (BCV)
Purity	≥85% (SDS-PAGE)
Sequence	MASLSGPISP INLEMFKPGV EELNPSKLLL LSNHQEGMLY PTILGSLELL SFKRERSLNL QRDKVCLLHQ ESQLLKL RGT GTDTT DVPLK QPMATSVNCC HDGIFTILEQ DRMPKTS MAP TLTESTGSLV TRLMSIPRLT FSIGTQVAMR LFRLGFRLAR YSLRV TILKA QEGLLLIPDL LHAHPVEPLV QDRAVEPILA IEPLPLV
Source	Baculovirus
Target Names	N
Protein Names	Recommended name: Protein I Alternative name(s): Accessory protein N2 N internal ORF protein Short name= IORF Protein in nucleocapsid ORF
Expression Region	1-207
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.