



Recombinant Mouse Histone-lysine N-methyltransferase SETD7 (Setd7)

Product Code	CSB-BP845030MO
Abbreviation	Setd7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VHL1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MDSDEEVVEE AVEGHLDDDG LPHGFCTVTY SSTDRFEGNF VHGEKNGRGK FFFFDGSTLE GYYVDDALQG QGVYTYEDGG VLQGTYVDGE LNGPAQEYDS DGRLIFKGQY KDNNRHGVCW IHYPDGGSLV GEVNEDGEMT GEKIAYVYPD QRTALYGKFI DGEMLEGKLA TLMATEEGRP HFEVTSGSSV YHFDKSTSSC ISSDALLPDP YESERVYVAD SLISSAGEGL FSKVAVGPNT VMSFYNGVRI THQEVDSRDW ALNGNTLSLD EETVIDVPEP YNHVSKEYCAS LGHKANHSFT PNCVYDLFVH PRFGPIKCIR TLRAVEAEEE LTVAYGYDHS PPGKSGPEAP EWYQVELKAF QATQQK
Source	Baculovirus
Target Names	Setd7
Protein Names	Recommended name: Histone-lysine N-methyltransferase SETD7 EC= 2.1.1.43 Alternative name(s): Histone H3-K4 methyltransferase SETD7 Short name= H3-K4-HMTase SETD7 SET domain-containing protein 7 SET7/9
Expression Region	1-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.