



Recombinant Rat Death domain-associated protein 6 (Daxx)

Product Code	CSB-MP845039RA
Abbreviation	Daxx
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VIB2
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MATDDSIIVL DDDDEDEAAA QPGPSNPASN PVSPGPEASG PSESHGDGGS SNSGSRKCYK LENEKLFEEF LELCKMQTSD HPEVVPFLHK LQQRQAQSAFL ASAEFRNLS RVLRSRNRNP AKLYVYINEL CTVLKAHSIK KKLNLAPAAS AESSGDNPPT DPPSDLTNTTE TTASEASRTR GSRRQIQRLQ LLLALYVAEI QRLQEKELDL SELDDPDSSY LQEARLKRKL IRLFGRLCDL KDCSSLTGKV IEQRIPYRGT RYPEVNRRIE RLINKPGPDT FPDYGDVLR VEGAATRHS GLPRQQLQLL AQDAFRDVG VRLQERRHLDL IYNFGCHLTD DYRPGVDPAL TDPTLARRLR ENRTLAMSRL DEVISKYAMM QDKSEEGERQ KRRARLLATS QSSDLPKASS DSGEGPSGVA SQEDPTTPKA ETEDEEDDEE SDDEEEEEEE EEEEATEDED EDLEQLQEDQ DDEEEEEEDN EDDKSPASPS PIFRRKEFSN PQKGGSPQEE QQERGLTGTP ASPLEASPGL PSTDAESSGE QLQERLLAGE SPVSQLSELD MEALPEETIP SPEERGISSS RRKSDSSLPT ILENGAAMVT STSFNGRVSS HPCRNASPPS KRFRKEKKQL GPGPLGNSYV KKQTMAQQDS GWKISVLSTP SSPLASVGPV ADSSTRVDSP SHELVTSSLC NPSPSLILQT PQSQSPRPCI YKTSVATQCD PEEIIVLSDS D
Source	Mammalian cell
Target Names	Daxx
Protein Names	Recommended name: Death domain-associated protein 6 Alternative name(s): Daxx
Expression Region	1-731
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a multifunctional protein that resides in multiple locations in the nucleus and in the cytoplasm. It interacts with a wide variety of proteins, such as apoptosis antigen Fas, centromere protein C, and transcription factor



erythroblastosis virus E26 oncogene homolog 1. In the nucleus, the encoded protein functions as a potent transcription repressor that binds to sumoylated transcription factors. Its repression can be relieved by the sequestration of this protein into promyelocytic leukemia nuclear bodies or nucleoli. This protein also associates with centromeres in G2 phase. In the cytoplasm, the encoded protein may function to regulate apoptosis. The subcellular localization and function of this protein are modulated by post-translational modifications, including sumoylation, phosphorylation and polyubiquitination. Alternative splicing results in multiple transcript variants.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.