



Recombinant Human B-cell CLL/lymphoma 7 protein family member C (BCL7C)

Product Code	CSB-YP845147HU
Abbreviation	BCL7C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WUZ0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAGRTVRAET RSRAKDDIKK VMATIEKVRR WEKRWVTVGD TSLRIFKWVP VVDPQEEERR RAGGGAERSR GRERRGRGAS PRGGGPLILL DLNDENSNQS FHSEGLQKG TEPSPGGTPQ PSRPVSPAGP PEGVPEEAQP PRLGQERDPG GITAGSTDEP PMLTKEEPVP ELLEAEAPEA YPVFEPVPPV PEAAQGDTEG SEGAPPLKRI CPNAPDP
Source	Yeast
Target Names	BCL7C
Protein Names	Recommended name: B-cell CLL/lymphoma 7 protein family member C
Expression Region	1-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is identified by the similarity of its product to the N-terminal region of BCL7A protein. The BCL7A protein is encoded by the gene known to be directly involved in a three-way gene translocation in a Burkitt lymphoma cell line. The function of this gene has not yet been determined.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.