



# Recombinant Human RNA polymerase-associated protein LEO1 (LEO1)

<b>Product Code</b>	CSB-EP845150HU-B
<b>Abbreviation</b>	LEO1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WVC0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ADMEDLFGS DADSEAERKD SDSGSDSDSD QENAASGSNA SGSESDQDER GDSGQPSNKE LFGDDSEDEG ASHHSGSDNH SERSDNRSEA SERSDHEDND PSDVDQHSGS EAPNDDEDEG HRSDGGSHHS EAEGSEKAHS DDEKWGREDK SDQSDDEKIQ NSDDEERAQG SDEDKLQNSD DDEKMQNTDD EERPQLSDDE RQQLSEEEKA NSDDERPVAS DNDDEKQNSD DEEQPQLSDE EKMQNSDDER PQASDEEHRH SDDEEEQDHK SESARGSDSE DEVLRMKRKN AIASDSEADS DTEVPKDNSG TMDLFGGADD ISSGSDGEDK PPTPGQPVDE NGLPQDQQEE EPIPETRIEV EIPKVNTDLG NDLYFVKLPN FLSVEPRPFD PQYYEDEFED EEMLDEEGRT RLKLVKVENTI RWRIRRDEEG NEIKESNARI VKWSDGMSL HLGNEVFDVY KAPLQGDHNN LFIRQGTGLQ GQAVFKTKLT FRPHSTDSAT HRKMTLSLAD RCSKTQKIRI LPMAGRDPEC QRTEMIKKEE ERLRASIRRE SQRRMREKQ HQRGLSASYL EPDRYDEEEE GEESISLAAI KNRYKGGIRE ERARIYSSDS DEGSEEDKAQ RLLKAKKLTS DEEGEPSGKR KAEDDDKANK KHKKYVISDE EEEDDD
<b>Source</b>	E.coli
<b>Target Names</b>	LEO1
<b>Protein Names</b>	Recommended name: RNA polymerase-associated protein LEO1 Alternative name(s): Replicative senescence down-regulated leo1-like protein
<b>Expression Region</b>	2-666
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.