



Recombinant Human Ankyrin repeat, SAM and basic leucine zipper domain-containing protein 1 (ASZ1)

Product Code	CSB-EP845159HU
Abbreviation	ASZ1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WWH4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAASALRGLP VAGGGESSES EDDGWEIGYL DRTSQKLKRL LPIEEKKEKF KKAMTIGDVS LVQELLDSGI SVDSNFQYGW TPLMYAASVA NAELVRVLLD RGANASFEKD KQSILITACS AHGSEEQILK CVELLSRNA DPNVACRRLM TPIMYAARDG HTQVVALLVA HGAEVNTQDE NGYTALTWAA RQGHKNIVLK LLELGANKML QTKDGKMPSE IAKRNKHHEI FNLLSFTLNP LEGKLQQLTK EDTICKILTT DSDREKDHIF SSYTAFGDLE VFLHGLGLEH MTDLLKERDI TLRHLLTMRE DEFTKNGITS KDQQKILAAL KELQVEEIQF GELSEETKLE ISGDEFNLNFL LKLNKQCGHL ITAVQNVITE LPVNSQKITL EWASPQNFTS VCEELVNNVE DLSEKVCKLK DLIQKLQNER ENDPHTIQLR EEVSTWNSRI LKRTAITICG FGFLLFICKL TFQRK
Source	E.coli
Target Names	ASZ1
Protein Names	Recommended name: Ankyrin repeat, SAM and basic leucine zipper domain-containing protein 1 Alternative name(s): Ankyrin-like protein 1 Germ cell-specific ankyrin, SAM and basic leucine zipper domain-containing protein
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.