



Recombinant Human Ankyrin repeat, SAM and basic leucine zipper domain-containing protein 1 (ASZ1)

| | |
|--------------------------|--|
| Product Code | CSB-MP845159HU |
| Abbreviation | ASZ1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8WWH4 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MAASALRGLP VAGGGESSES EDDGWEIGYL DRTSQKLRKRL LPIEEKKEKF KKAMTIGDVS LVQELLDSCI SVDSNFQYGW TPLMYAASVA NAELVRVLLD RGANASFEKD KQSILITACS AHGSEEQILK CVELLSRNA DPNVACRRLM TPIMYAARDG HTQVVALLVA HGAEVNTQDE NGYTALTWAA RQGHKNIVLK LLELGANKML QTKDGMKMPSE IAKRNRKHEI FNLLSFTLNP LEGKLQQLTK EDTICKILTT DSDREKDHIF SSYTAFGDLE VFLHGLGLEH MTDLLKERDI TLRHLLTMRE DEFTKNGITS KDQQKILAL KELQVEEIQF GELSEETKLE ISGDEFNLFL LKLNKQCGHL ITAVQNVITE LPVNSQKITL EWASPQNFTS VCEELVNNVE DLSEKVCKLK DLIQKLQNER ENDPHTIQLR EEVSTWNSRI LKRITAITICG FGFLLFICKL TFQRK |
| Source | Mammalian cell |
| Target Names | ASZ1 |
| Protein Names | Recommended name: Ankyrin repeat, SAM and basic leucine zipper domain-containing protein 1 Alternative name(s): Ankyrin-like protein 1 Germ cell-specific ankyrin, SAM and basic leucine zipper domain-containing protein |
| Expression Region | 1-475 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.