



Recombinant Human Secretogranin-3 (SCG3)

Product Code	CSB-BP845166HU
Abbreviation	SCG3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WXD2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	F PKPGGSQDKS LHNRELSAER PLNEQIAEAE EDKIKKTYPP ENKPGQSNYS FVDNLNLLKA ITEKEKIEKE RQSIRSSPLD NKLNVEDVDS TKNRKLIDY DSTKSGLDHK FQDDPDGLHQ LDGTPLEAED IVHKIAARIY EENDRAVFDK IVSKLLNLGL ITESQAHTLE DEVAEVLQKL ISKEANNYEE DPNKPTSWTE NQAGKIPEKV TPMAAIQDGL AKGENDETVS NTLTLNGLE RRTKTYSEDN FEELQYFPNF YALLKSIDSE KEAKEKETLI TIMKTLIDFV KMMVKYGTIS PEEGVSYLEN LDEMIALQTK NKLEKNATDN ISKLFAPSE KSHEETDSTK EEAAKMEKEY GSLKDSTKDD NSNPGGKTDE PKGKTEAYLE AIRKNIEWLK KHDKKGNKED YDLSKMRDFI NKQADAYVEK GILDKEEAEA IKRIYSSL
Source	Baculovirus
Target Names	SCG3
Protein Names	Recommended name: Secretogranin-3 Alternative name(s): Secretogranin III Short name= SgIII
Expression Region	20-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a member of the chromogranin/secretogranin family of neuroendocrine secretory proteins. Granins may serve as precursors for biologically active peptides. Some granins have been shown to function as helper proteins in sorting and proteolytic processing of prohormones; however, the function of this protein is unknown. Two transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.