



Recombinant Human Microprocessor complex subunit DGCR8 (DGCR8)

Product Code	CSB-BP845175HU
Abbreviation	DGCR8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WYQ5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	METDESPSPLPCGPAGEAVMESRARPFQALPREQSPPPPLQTSSGAEVMDV GSGGDGQSE LPAEDPFNFYGASLLSKGSFSKGRLLIDPNCSGHSPRTARHAPAVRKFSPDLK LLKDVKI SVSFTESCRSKDRKVLYTGAERDVRAECGLLLSPVSGDVHACPFGGSVGDGV GIGGESAD KKDEENELDQEKRVEYAVLDELEDFTDNLELDEEGAGGFTAKAIVQRDRVDEE ALNFPYE DDFDNDVDALLEEGLCAPKKRRTEEKYGGSDHPSDGETSVQPMMTKIKTVL KSRGRPPT EPVL
Source	Baculovirus
Target Names	DGCR8
Protein Names	Recommended name: Microprocessor complex subunit DGCR8 Alternative name(s): DiGeorge syndrome critical region 8
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein(Isoform 2)
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.