



Recombinant Human Kelch-like protein 6 (KLHL6)

Product Code	CSB-BP845176HU
Abbreviation	KLHL6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WZ60
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLMAGQRGAW TMGDVVEKSL EGPLAPSTDE PSQKTGDLVE ILNGEKVKFD DAGLSLILQN GLETLRMENA LTDVILCVDI QEFSCHRVVL AAASNYFRAM FCNDLKEKYE KRIIIKGVDA ETMHTLLDYT YTSKALITKQ NVQRVLEAAN LFQFLRMVDA CASFLTEALN PENCVGILRL ADTHSLDSLK KQVQSYIIQN FVQILNSEEF LDLPVDTLHH ILKSDDLYVT EEAQVFETVM SWVRHKPSER LCLLPYVLEN VRLPLDPWY FVETVEADPL IRQCPEVFPL LQEARMYHLS GNEIISERTK PRMHEFQSEV FMIIGGCTKD ERFVAEVTCL DPLRRSRLEV AKLPLTEHEL ESENKKWVEF ACVTLKNEVY ISGGKETQHD VWKYNSSINK WIQIEYLNIG RWRHKMVVLG GKVYVIGGFD GLQRINNVET YDPFHNCWSE AAPLLVHVSS FAATSHKKKL YVIGGGPNGK LATDKTQCYD PSTNKWLSKA AMPVEAKCIN AVSFRDRIYV VGGAMRALYA YSPLEDSWCL VTQLSHERAS CGIAPCNNRL YITGGRDEKN EVIATVLCWD PEAQKLTEEC VLPRGVSHHG SVTIRKSYTH IRRIVPGAVS V
Source	Baculovirus
Target Names	KLHL6
Protein Names	Recommended name: Kelch-like protein 6
Expression Region	1-621
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.