



# Recombinant Human Kelch-like protein 6 (KLHL6)

<b>Product Code</b>	CSB-YP845176HU
<b>Abbreviation</b>	KLHL6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WZ60
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLMAGQRGAW TMGDVVEKSL EGPLAPSTDE PSQKTGDLVE ILNGEKVKFD DAGLSLILQN GLETLRMENA LTDVILCVDI QEFSCHRVVL AAASNYFRAM FCNDLKEKYE KRIIIKGVDA ETMHTLLDYT YTSKALITKQ NVQRVLEAAN LFQFLRMVDA CASFLTEALN PENCVGILRL ADTHSLDSLK KQVQSYIIQN FVQILNSEEF LDLPVDTLHH ILKSDDLYVT EEAQVFETVM SWVRHKPSE LCLLPYVLEN VRLPLDPWY FVETVEADPL IRQCPEVFP LQEARMYHLS GNEIISERTK PRMHEFQSEV FMIIGGCTKD ERFVAEVTCL DPLRRSRLEV AKLPLTEHEL ESENKKWVEF ACVTLKNEVY ISGGKETQHD VWKYNSSINK WIQIEYLNIG RWRHKMVVLG GKVVYVIGGFD GLQRINNVET YDPFHNCWSE AAPLLVHVSS FAATSHKKKL YVIGGGPNGK LATDKTQCYD PSTNKWLSKA AMPVEAKCIN AVSFRDRIYV VGGAMRALYA YSPLEDSWCL VTQLSHERAS CGIAPCNNRL YITGGRDEKN EVIATVLCWD PEAQKLTEEC VLPRGVSHHG SVTIRKSYTH IRRIVPGAVS V
<b>Source</b>	Yeast
<b>Target Names</b>	KLHL6
<b>Protein Names</b>	Recommended name: Kelch-like protein 6
<b>Expression Region</b>	1-621
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.