



Recombinant Putative xanthine dehydrogenase YagS FAD-binding subunit (yagS)

Product Code	CSB-MP845203EOD
Abbreviation	yagS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8X6J0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MKAFTYERVN TPAEAALSAQ RVPGAKFIAG GTNLLDLMKL EIETPTHLID VNGLGLDKIE VTDAGGLRIG ALVRNTDLVA HERVRRDYAV LSRALLAGAS GQLRNQATTA GNLLQRTRCP YFYDTNQCPCN KRLPGSGCAA LEGFSRQHAV VGVSEACIAT HPSDMAVAMR LLDAVVETIT PEGKTRSITL ADFYHPPGKT PHIETALLPG ELIVAVTLPP PLGGKHIYRK VRDRASYTFA LVSVAIIQP DGSGRVALGG VAHKPWRIEA ADAQLSQGAQ AVYDALFASA HPTAENTFKL LLAKRTLASV LAEARAQA
Source	Mammalian cell
Target Names	yagS
Protein Names	Recommended name: Putative xanthine dehydrogenase YagS FAD-binding subunit EC= 1.17.1.4
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.