



Recombinant Beet curly top virus Protein C2 (C2, L2)

Product Code	CSB-BP845964BFX
Abbreviation	C2, L2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91J26
Product Type	Recombinant Protein
Immunogen Species	Beet curly top virus (strain California/Logan) (BCTV)
Purity	>85% (SDS-PAGE)
Sequence	MENHVSLKVV SPALYYAIQD LRAHTNNFLK NQKMKPLSPG HYIIQPSANS KVRSLITKQQ HPRKVTLPEN CHFTIHHECN RGFSHRGTYT SPSGKFKRGI RECTESTVYE TPMVREIRAN LSTEDTNPIQ LQPPESSVESS QVLDRANDNR IEQDIDWTPF LEGLEKETRD ILG
Source	Baculovirus
Target Names	C2
Protein Names	Recommended name: Protein C2 Alternative name(s): Protein L2
Expression Region	1-173
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.