



Recombinant Mouse Interferon-inducible protein AIM2 (Aim2)

| | |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-YP845977MO |
| Abbreviation | Aim2 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q91VJ1 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MESEYREMLL LTGLDHITTE ELKRFKYFAL TEFQIARSTL DVADRTELAD HLIQSAGAAS AVTKAINIFQ KLNMYMHIANA LEEKKKEAER KLMTNTKKRG TQKVENRSQA ENCSAASATR SDNDFKEQAA TEVCPQAKPQ KKQMVAEQEA IREDLQKDPL VVTVLKAINP FECETQEGRQ EIFHATVATE TDFFFVKVLN AQFKDKFIPK RTIKISNYLW HSNFMEVTSS SVVVDVESNH EVPNNVVKRA RETPRISKLK IQPCGTIVNG LFKVQKITEE KDRVLYGIHD KTGTMEVLVL GNPSKTKCEE GDKIRLTFE VSKNGVKIQL KSGPCSFFKV IKAAKPKTDM KSVE |
| Source | Yeast |
| Target Names | Aim2 |
| Protein Names | Recommended name: Interferon-inducible protein AIM2 Alternative name(s): Interferon-inducible protein 210 Short name= lfi-210 Interferon-inducible protein p210 |
| Expression Region | 1-354 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Target Details | AIM2 is a member of the IFI20X /IFI16 family. It plays a putative role in tumorigenic reversion and may control cell proliferation. Interferon-gamma induces expression of AIM2. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.