



Recombinant Rat F-box only protein 6 (Fbxo6)

Product Code	CSB-YP846061RA
Abbreviation	Fbxo6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q923V4
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MVNINELPEN ILLELFTHVP APQLLRNCRL VCSSLWRDLID VMTLWKRKSL REGFVTKDRD EPVDDWKIFY ILCSLQRNLL RNPCAENLR SWRIDSNNGD EWKVESLPGD HGTSFPDTKV KKYFVTSYGM CLKSQMVDLK AEGYSEKLLD TVRPDIVVKD WFAPRADCGC TYHLRVQLAS ADYIVLASFE PPPVTIEQWN DASWQEISHT FSNYPPGVRH ILFQHGGKDT QFWKGWYGPR VTNSSIIVSH RTAKNPAPAR TLPEEDTSNR RKILSFGSWE DLSP
Source	Yeast
Target Names	Fbxo6
Protein Names	Recommended name: F-box only protein 6 Alternative name(s): F-box only protein 6b F-box protein that recognizes sugar chains 2 F-box/G-domain protein 2
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. This protein belongs to the Fbxs class, and its C-terminal region is highly similar to that of rat NFB42 (neural F Box 42 kDa) which may be involved in the control of the cell cycle.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.