



# Recombinant *Cavia porcellus* Interleukin-12 subunit beta (IL12B)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP846073GU  |
| <b>Abbreviation</b>      | IL12B   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q924V5  |
| <b>Storage Buffer</b>    | Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0  |
| <b>Product Type</b>      | Recombinant Proteins  |
| <b>Immunogen Species</b> | <i>Cavia porcellus</i> (Guinea pig)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MWELKKDV YVVELDWHTD APGETVVLTC NTAEDGITW TSDRKSDILG<br>SGKTLTIQVK EFEDAGGYTC HKGGEVLSRS QLLLHKKEDE IWSTDILKEQ<br>KGSNGKTLK CEARSYSGRF TCWWLTAFGT DVKFSVKGSR GSSDPSGVTC<br>GEAERVSGDN QEYKYSVECQ EDSACPTAEE SLPIEVVDA IHKFKYENYT<br>SSFYIRDIK PDPPKNLQLK PSVNSQQVEV SWEYPDTWST PHSYFSLTFL<br>VQTHGKNKNR REKKYELFTD KTSATVSCHK ISKVEVRARD RYYSSSWSEW<br>ASVSCSEVSV SR   |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | IL12B   |
| <b>Protein Names</b>     | Recommended name: Interleukin-12 subunit beta Short name=<br>IL-12B Alternative name(s): Cytotoxic lymphocyte maturation factor 40 kDa<br>subunit Short name= CLMF p40 IL-12 subunit p40  |
| <b>Expression Region</b> | 23-332  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.