



# Recombinant Human Transcription factor GATA-6 (GATA6)

<b>Product Code</b>	CSB-BP846113HU
<b>Abbreviation</b>	GATA6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q92908
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MALTDGGWCL PKRFGAAGAD ASDSRAFFAR EPSTPPSPIS SSSSSCSRGG ERGGGASNC GTPQLDTEAA AGPPARLLLL SSYASHPFGA PHGPSAPGVA GPGGNLSSWE DLLLFTDLQ AATASKLLWS SRGAKLSPFA PEQPEEMYQT LAALSSQGPA AYDGAPGGFV HSAAAAAAAAA AAASSPVYVP TTRVGSMLPG LPYHLQGSGS GPANHAGGAG AHPGWPQASA DSPPYGSGGG AAGGGAAGPG GAGSAAAHVS ARFPYSPSP MANGAAREPG GYAAAGSGGA GGVSGGGSSL AAMGGREPQY SSLSAARPLN GTYHHHHHHH HHHPSYSPY VGAPLTPAWP AGPFETPV LH SLQSRAGAPL PVPRGPSADL LEDLSESREC VNCGSIQTPL WRRDGTGHYL CNACGLYSKM NGLSRPLIKP QKRVPSRRL GLSCANCHTT TTTLWRRNAE GEPVCNACGL YMKLHGVRP LAMKKEGIQT RKRKPKNINK SKTCSGNSNN SIPMTPTSTS SNSDDCSKNT SPTTQPTASG AGAPVMTGAG ESTNPENSEL KYSGQDGLYI GVSLASPAEV TSSVRPDSWC ALALA
<b>Source</b>	Baculovirus
<b>Target Names</b>	GATA6
<b>Protein Names</b>	Recommended name: Transcription factor GATA-6 Alternative name(s): GATA-binding factor 6
<b>Expression Region</b>	1-595
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.