



Recombinant Human Transcription factor GATA-6 (GATA6)

Product Code	CSB-EP846113HU
Abbreviation	GATA6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92908
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MALTDGGWCL PKRFGAAGAD ASDSRAFFAR EPSTPPSPIS SSSSSCSRGG ERPGGGASNC GTPQLDTEAA AGPPARSLLL SSYASHPFGA PHGPSAPGVA GPGGNLSSWE DLLLFTDLQ AATASKLLWS SRGAKLSPFA PEQPEEMYQT LAALSSQGPA AYDGAPGGFV HSAAAAAAAAA AAASSPVYVP TTRVGSMLPG LPYHLQGSGS GPANHAGGAG AHPGWPQASA DSPPYGSGGG AAGGGAAGPG GAGSAAAHVS ARFPYSPSP MANGAAREPG GYAAAGSGGA GGVSGGGSSL AAMGGREPQY SSLSAARPLN GTYHHHHHHH HHHPSYSPY VGAPLTPAWP AGPFETPVLH SLQSRAGAPL PVPRGPSADL LEDLSESREC VNCGSIQTPL WRRDGTGHYL CNACGLYSKM NGLSRPLIKP QKRVPSSRRL GLSCANCHTT TTTLWRRNAE GEPVCNACGL YMKLHGVRP LAMKKEGIQT RKRKPKNINK SKTCGNSNN SIPMTPST SNSDDCSKNT SPTTQPTASG AGAPVMTGAG ESTNPENSEL KYSGQDGLYI GVSLASPAEV TSSVRPDSWC ALALA
Source	E.coli
Target Names	GATA6
Protein Names	Recommended name: Transcription factor GATA-6 Alternative name(s): GATA-binding factor 6
Expression Region	1-595
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.