



Recombinant *Synechococcus elongatus* Glycogen synthase (glgA)

Product Code	CSB-EP846311FPY-B
Abbreviation	glgA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q935Y7
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MRILFVAAEC APFAKVGGMG DVVGSPLPKVL KALGH DVRIF MPYYGFLNSK LDIPAEPIWW GYAMFNHFAV YETQLPGSDV PLYLMGHPAF DPHRIYSGED EDWRFTFFAN GAAEFSWNYW KPQVIHCHDW HTGMIPVWMH QSPDISTVFT IHNLAYQGPW RWKLEKITWC PWYMQGDSTM AAALLYADRV NTVSPTYAQQ IQTPTYGEKL EGLLSFISGK LSGILNGIDV DSYNPA TDTR IVANYDRDTL DKRLNNKLAL QKEMGLEVNP DRFLIGFVAR LVEQKGIDLL LQILDRFLSY SDAQFVVLGT GERYYETQLW ELATRYPGRM STYLMYDEGL SRRIYAGSDA FLVPSRFEPG GITQMLALRY GSVPIVRRTG GLVDTV FHHD PRHAEGNGYC FDRYEPLDLY TCLVRAWESY QYQPQWQKLQ QRGMAVDLSW KQSAIAYEQL YAEAIGLPID VLQEA
Source	E.coli
Target Names	glgA
Protein Names	Recommended name: Glycogen synthase EC= 2.4.1.21 Alternative name(s): Starch [bacterial glycogen] synthase
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.