



# Recombinant Arabidopsis thaliana Ribose-phosphate pyrophosphokinase 3, chloroplastic (PRS3)

<b>Product Code</b>	CSB-BP846364DOA
<b>Abbreviation</b>	PRS3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q93Z66
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	T ASVSNRCIKC GVRSLNHSG HRSLDFLSNG DPISLINPNS SSPITMAAAT SESGSKSSKR VCLFHSDETR DLAERIVAKS DCIELRSINW KKFDDGFPNL FIQNAQGIRG QHVAFLASFS SPAVIFEQLS VIYALPKLFV SSFTLVLPFF PTGTSERMED EGDVATAFTL ARILSNIPTS RGGPTSLVTF DIHALQERFY FGDTILPCFE SGIPLLKSR LQSLPDSNIS IAFPDDGAWK RFHKQLQHYP TIVCNKVRMG DKRIVRIKEG DAEGRHVVIV DDLVQSGGTL IECQKVLAAH GAAKISAYVT HGIFPRSSWK RFKLDTKGDP AEGLSYFWIT DSCGMTVKEV MNKPPFEVLS LAGSIASALQ V
<b>Source</b>	Baculovirus
<b>Target Names</b>	PRS3
<b>Protein Names</b>	Recommended name: Ribose-phosphate pyrophosphokinase 3, chloroplastic EC= 2.7.6.1 Alternative name(s): Phosphoribosyl pyrophosphate synthase 3
<b>Expression Region</b>	40-411
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.