



# Recombinant Arabidopsis thaliana Tubby-like F-box protein 5 (TULP5)

<b>Product Code</b>	CSB-EP846386DOA-B
<b>Abbreviation</b>	TULP5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q944S3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSFLSIVRDV RDTVGSFSRR SFDVRVSNGT THQRSKSHGV EAHIEDLIVI KNTRWANLPA ALLRDVMKKL DESESTWPAR KQVVACAGVC KTWRLMCKDI VKSPEFSGKL TFPVSLKQPG PRDGIIQCYI KRDKSNMTYH LYLSLSPAIL VESGKFLLSA KRSRRATYTE YVISMDADNI SRSSSTYIGK LKSNFLGTKF IVYDTAPAYN SSQILSPPNR SRSFNSKKVS PKVPSGSYNI AQVTYELNLL GTRGPRRMNC IMHSIPSLAL EPGGTVPSPQ EFLQRSLDES FRSIGSSKIV NHSGDFTRPK EEEGKVRPLV LKTKPPRWLQ PLRCWCLNFK GRVTVASVKN FQLMSAATVQ PGSGSDGGAL ATRPSLSPQQ PEQSNHDKII LHFGKVGKDM FTMDYRYPLS AFQAF AISLS TFDTKLACE
<b>Source</b>	E.coli
<b>Target Names</b>	TULP5
<b>Protein Names</b>	Recommended name: Tubby-like F-box protein 5 Short name= AtTLP5
<b>Expression Region</b>	1-429
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.