



Recombinant Arabidopsis thaliana Tubby-like F-box protein 5 (TULP5)

Product Code	CSB-EP846386DOA
Abbreviation	TULP5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q944S3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSFLSIVRDV RDTVGSFSRR SFDVRVSNGT THQRSKSHGV EAHIEDLIVI KNTRWANLPA ALLRDVMKKL DESESTWPAR KQVVACAGVC KTWRLMCKDI VKSPEFSGKL TFPVSLKQPG PRDGIIQCYI KRDKSNMTYH LYLSLSPAIL VESGKFLLSA KRSRRATYTE YVISMDADNI SRSSSTYIGK LKSNFLGTKF IVYDTAPAYN SSQILSPPNR SRSFNSKKVS PKVPSGSYNI AQVTYELNLL GTRGPRRMNC IMHSIPSLAL EPGGTVPSPQ EFLQRSLDES FRSIGSSKIV NHSGDFTRPK EEEGKVRPLV LKTKPPRWLQ PLRCWCLNFK GRVTVASVKN FQLMSAATVQ PGSGSDGGAL ATRPSLSPQQ PEQSNHDKII LHFGKVGKDM FTMDYRYPLS AFQAFASLS TFDTKLACE
Source	E.coli
Target Names	TULP5
Protein Names	Recommended name: Tubby-like F-box protein 5 Short name= AtTLP5
Expression Region	1-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.