



Recombinant *Oryza sativa* subsp. japonica Probable E3 ubiquitin-protein ligase XBOS31 (XBOS31)

Product Code	CSB-EP846418OFG-B
Abbreviation	XBOS31
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94CT7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGHGLSCSRD TDEYDLFRAA QLGDHALSA LLAADPALAR RATVYDRFTA LHIAANGRL QVLSMLLDRD GDVDVLSRKK QTPLMVAAMR GNTECVVRL RGGANVLTFD SPRARTCLHH AAYYGHAEC QAILGAAAQA QGPVAASWGF ARFVNVRDER GATPLHLAAR HARASCVRL LDKGAIVSAP TAVYGFPGST ALHLAARAGS MECIRELLAW GADRLQRDSA GRIAYAVAMR RGHRACAALL NPAAAEPIVW PSPLKFIGEL EADAKALLEA ALMEANRERE KRILHGSDIN IKGGDEEEES EDEEEACNIC FEQACSM EVK ECGHQMCAAC TLAICCHSKP NPKTLLLHPP ACPFCRTTIS RLVVATTNSN KTNSRRRSRS RSSSFKGGLS SAMGSFSRIG RGSRLVVDG SSVGELADKP DHDFSSVAAA AAICDT
Source	E.coli
Target Names	XBOS31
Protein Names	Recommended name: Probable E3 ubiquitin-protein ligase XBOS31 EC= 6.3.2.- Alternative name(s): Ankyrin repeat domain and RING finger-containing protein XBOS31 XB3 protein homolog 1
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.