



# Recombinant Bovine Pigment epithelium-derived factor (SERPINF1)

<b>Product Code</b>	CSB-YP846450BO
<b>Abbreviation</b>	SERPINF1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q95121
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	NAGQEAGSLT PESTGAPVEE EDPFFKVPVN KLAAAVSNFG YDLRVRSGE SPTANVLLSP LSVATALSAL SLGAEQRTES NIHRALYYDL ISNPDIHGTY KDLLASVTAP QKNLKSASRI IFERKLRIKA SFIPPLEKSY GTRPRILTGN SRVDLQEINN WVQAQMKGV ARSTREMPSE ISIFLLGVAY FKGQWVTKFD SRKTSLEDFY LDEERTVKVP MMSDPQAVLR YGLDSDLNCK IAQLPLTGST SIIFFLPQKV TQNLTLIEES LTSEFIHDID RELKTVQAVL TIPKLKLSYE GELTKSVQEL KLQSLFDAPD FSKITGKPIK LTQVEHRVGF EWNEDGAGTN SSPGVQPARL TFPLDYHLNQ PFIFVLRDTD TGALLFIGKI LDPRGT
<b>Source</b>	Yeast
<b>Target Names</b>	SERPINF1
<b>Protein Names</b>	Recommended name: Pigment epithelium-derived factor Short name= PEDF Alternative name(s): Serpin F1
<b>Expression Region</b>	21-416
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.