



Recombinant *Macaca fascicularis* Serine/threonine-protein phosphatase 2A 55 kDa regulatory subunit B gamma isoform (PPP2R2C)

Product Code	CSB-EP846506MOV-B
Abbreviation	PPP2R2C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q95LP0
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Purity	>85% (SDS-PAGE)
Sequence	MGLSFFSKHL PIQEGQPWAS KTPADIISTV EFNHTGELLA TGDKGGRVVI FQREPESKNA PHSQGEYDVY STFQSHEPEF DYLSKLEIEE KINKIKWLPQ QNAAHSLST NDKTIKLWKI TERDKRPEGY NLKDEEGKLG DLSTVTSLQV PVLKPMMLMV EVSPRRIFAN GHTYHINSIS VNSDCETYMS ADDLRINLWH LAITDRSFNI VDIKPANMED LTEVITASEF HPHHCNLFVY SSSKGSLRLC DMRAAALCDK HSKLFEEPED PSNRSFFSEI ISSVSDVKFS HSGRYMLTRD YLTVKVWDLN MEARPIETYQ VHDYLRSLKC SLYENDCIFD KFECAWNGSD SVIMTGAYNN FFRMFDRNTK RDVTLEASRE SSKPRAVLKP RRVCVGGKRR RDDISVDSL DFTKKILHTAW HPAENIIAIA ATNNLYIFQD KVNSDMH
Source	<i>E.coli</i>
Target Names	PPP2R2C
Protein Names	Recommended name: Serine/threonine-protein phosphatase 2A 55 kDa regulatory subunit B gamma isoform Alternative name(s): PP2A subunit B isoform B55-gamma PP2A subunit B isoform PR55-gamma PP2A subunit B isoform R2-gamma PP2A subunit B
Expression Region	1-447
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.