



Recombinant *Drosophila melanogaster* Attacin-C (AttC)

Product Code	CSB-EP846516DLU
Abbreviation	AttC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q95NH6
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	QRPYTQP LIYYPPPPTP PRIYRARRQV LGGSLTSNPS GGADARLDLS KAVGTPDHHV IGQVFAAGNT QTKPVSTPVT SGATLGYNNH GHGLELTKTH TPGVRDSFQQ TATANLFNNG VHNLDAKFA SQNQLANGFK FDRNGAALDY SHIKGHGATL THANIPGLGK QLELGGRANL WQSQDRNTRL DLGSTASKWT SGPFKGQTDL GANLGLSHYF G
Source	E.coli
Target Names	AttC
Protein Names	Recommended name: Attacin-C Cleaved into the following chain: 1. Immune-induced peptide 16 Short name= 2. DIM-16 Alternative name(s): MPAC
Expression Region	24-241
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.