



Recombinant Beta-lactamase-like protein 2 homolog (Y53F4B.39)

Product Code	CSB-MP846520CXY
Abbreviation	Y53F4B.39
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q95Q18
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MPSLTHVEPI EKLSDAVTRI LGHNPGPFTL QGTNTYLLGT GAKKILVDTG EPNVTEYISA LKSVLASTNS HIEYIVITHW HGDHVGIDN ITDEILDKKK IPIYKMKRDK DEGVERFHVY DDGFEVAVDG ATLKLIATPG HTADHFSLWL QEERALFSGD CILGEGTTVF EDLHDYMTSL QKIKDLNATR IYPGHGPVID KVVEKVDEYI EHRMKREREI IKVLKEHEEI TSMDVTNQQVY ADSPWAVRLA ALNNVKLVLK KLCKDGVVEN PHFETFKWIG GSSSGEEKKN ESSNL
Source	Mammalian cell
Target Names	Y53F4B.39
Protein Names	Recommended name: Beta-lactamase-like protein 2 homolog EC= 3.-.-.-
Expression Region	1-295
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.