



# Recombinant Human UPF0498 protein KIAA1191 (KIAA1191)

<b>Product Code</b>	CSB-BP846582HU
<b>Abbreviation</b>	KIAA1191
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96A73
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASRQPEVPA LEASAPLGKM SLPIGIYRRA VSYDDTLEDP APMTPPPSDM GSVPWKPVIP ERKYQHLAKV EEGEASLPSP AMTLSSAIDS VDKVPVVKAK ATHVIMNSLI TKQTQESIQH FERQAGLRDA GYTPHKGLTT EETKYLRVAE ALHKLKLQSG EVTKEERQPA SAQSTPSTTP HSSPKQRPRG WFTSGSSTAL PGPNPSTMDS GSGDKDRNLS DKWSLFGPRS LQKYDSGSFA TQAYRGAQKP SPLELIRAQA NRMAEDPAAL KPPKMDIPVM EGKKQPPRAH NLKPRDLNVL TPTGF
<b>Source</b>	Baculovirus
<b>Target Names</b>	KIAA1191
<b>Protein Names</b>	Recommended name: UPF0498 protein KIAA1191 Alternative name(s): Brain-derived rescue factor p60MONOX
<b>Expression Region</b>	1-305
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.