



Recombinant Human F-box/LRR-repeat protein 8 (FBXL8)

Product Code	CSB-MP846597HU
Abbreviation	FBXL8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96CD0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAEPGEGLP EVLALIFRHL SLRDRAAAAR VCRAWAAAT CSAVWHDTKI SCECELEGML PPYLSACLDH IHNLRLEFEP SRKPSRRAAI ELLMVLAGRA PGLRGLRLEC RGEKPLFDAG RDVLEAVHAV CGAASQLRHL DLRRLSFTLD DALVLQAARS CPELHSLFLD NSTLVGSGVP GSVLELLEAC PRLRALGLHL ASLSHAILEA LAAPDRAPFA LLALRCACPE DARASPLPNE AWVALRRRHP GLAVELELEP ALPAESVTRV LQPAVPVAAL RLNLSGDTVG PVRFAAHHYA ATLCALEVRA AASAEINAAL EELAARCAAL REVHCFCVVS HSVLDAFRAH CPRLRTYTLK LTREPHPWRP TLVA
Source	Mammalian cell
Target Names	FBXL8
Protein Names	Recommended name: F-box/LRR-repeat protein 8 Alternative name(s): F-box and leucine-rich repeat protein 8 F-box protein FBL8
Expression Region	1-374
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. This protein belongs to the Fbls class. It shares 78% sequence identity with the mouse protein.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.