



Recombinant Human Probable proline racemase (C14orf149)

Product Code	CSB-BP846613HU
Abbreviation	C14orf149
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96EM0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MESALAVPRL PPHDPGTPVL SVVDMHTGGE PLRIVLAGCP EVSGPTLLAK RRYMRQHLDH VRRRLMFEPGR GHRDMYGAVL VPSELPDAHL GVLFLHNEGY SSMCGHAVLA LGRFALDFGL VPAPPAGTRE ARVNIHCPCG LVTAFCVACED GRSHGPVRFH SVPFVFLATD LMVDVPGHGK VMVDIAYGGA FYAFVTAEKL GLDICS AKTR DLVDAASAVT EAVKAQFKIN HPDSEDLAFL YGTILT DGKD AYTKEPTTNI CVFADEQVDR SPTGSGVTAR IALQYHKGLL ELNQMRAFKS SATGSVFTGK AVREAKCGDF KAVIVEVSGQ AHYTG TASFI IEDDDPLRDG FLLK
Source	Baculovirus
Target Names	L3HYPDH
Protein Names	Recommended name: Probable proline racemase EC= 5.1.1.4
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.