



Recombinant Human Multivesicular body subunit 12A (FAM125A)

Product Code	CSB-MP846615HU
Abbreviation	FAM125A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96EY5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDPVPGTDSA PLAGLAWSSA SAPPPRGFSA ISCTVEGAPA SFGKSFAQKS GYFLCLSSLG SLENPQENVV ADIQIVVDKS PLPLGFSPVC DPMDSKASVS KKKRMCVKLL PLGATDTAVF DVRLSGKTKT VPGYLRIGDM GGFAIWCKKA KAPRPVPKPR GLSRDMQGLS LDAASQPSKG GLLERTASRL GSRASTLRRN DSIYEASSLY GISAMDGVPF TLHPRFEGKS CSPLAFSAFG DLTIKSLADI EEEYNYGFVV EKTAARLPP SVS
Source	Mammalian cell
Target Names	MVB12A
Protein Names	Recommended name: Multivesicular body subunit 12A Alternative name(s): CIN85/CD2AP family-binding protein ESCRT-I complex subunit MVB12A Protein FAM125A
Expression Region	1-273
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.