



Recombinant Human Phytanoyl-CoA hydroxylase-interacting protein-like (PHYHIPL)

Product Code	CSB-EP846618HU-B
Abbreviation	PHYHIPL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96FC7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEVPRLDHAL NSPTSPCEEV IKNLSLEAIQ LCDRDGNKSQ DSGIAEMEEL PVPHNIKISN ITCDSFKISW EMDSKSKDRI THYFIDLNKK ENKNSNKFKH KDVPTKLVAK AVPLPMTVRG HWFLSPRTEY TVAVQTASKQ VGDYVVSEW SEIIEFCTAD YSKVHLTQLL EKAEVIAGRM LKFSVFYRNQ HKEYFDYVRE HHGNAMQPSV KDNSGSHGSP ISGKLEGIFF SCSTEFNTGK PPQDSPYGRY RFEIAAEKLF NPNTNLYFGD FYCMYTAYHY VILVIAPVGS PGDEFCKQRL PQLNSKDNKF LTCTEEDGVL VYHHAQDVIL EVIYTPVDL SVGTVAEITG HQLMSLSTAN AKKDPSCKTC NISVGR
Source	E.coli
Target Names	PHYHIPL
Protein Names	Recommended name: Phytanoyl-CoA hydroxylase-interacting protein-like
Expression Region	1-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.