



Recombinant Human Charged multivesicular body protein 6 (CHMP6)

Product Code	CSB-MP846622HU
Abbreviation	CHMP6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96FZ7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GNLFGRRKKQ SRVTEQDKAI LQLKQQRDKL RQYQKRIAQQ LERERALARQ LLRDGRKERA KLLLKKKRYQ EQLLDRTENQ ISSLEAMVQS IEFTQIEMKV MEGLQFGNEC LNKMHQVMSI EEVERILDET QEAVEYQRQI DELLAGSFTQ EDEDAILLEEL SAITQEQIEL PEVPSEPLPE KIPENVPVKA RPRQAELVAA S
Source	Mammalian cell
Target Names	CHMP6
Protein Names	Recommended name: Charged multivesicular body protein 6 Alternative name(s): Chromatin-modifying protein 6 Vacuolar protein sorting-associated protein 20 Short name= Vps20 Short name= hVps20
Expression Region	2-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.