



Recombinant Human Lethal (3)malignant brain tumor-like protein 3 (L3MBTL3)

Product Code	CSB-BP846646HU
Abbreviation	L3MBTL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96JM7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTESASSTSG QEFDFVFSVMD WKDGVGTLPG SDLKFRVNEF GALEVITDEN EMENVKATA TTTWMVPTAQ EAPTSPSSR PVFPPAYWTS PPGCPTVFSE KTGMPFRLKD PVKVEGLQFC ENCCQYGNVD ECLSGGNYCS QNCARHIKDK DQKEERDVEE DNEEEDPKCS RKKKPKLSLK ADTKEDGEER DDEMENKQDV RILRGSQRAR RKRRGDSAVL KQGLPPKGKK AWCWASYLEE EKAVAVPAKL FKEHQSFYIN KNGFKVGMKL EGVDPHQSV YCVLTVAEVC GYRIKLFHFDG YSDCYDFWVN ADALDIHPVG WCEKTGHKLH PPKGYKEEEF NWQTYLKTCK AQAAPKSLFE NQNITVIPSG FRVGMKLEAV DKKNPSFICV ATVTDMVDNR FLVHFDNWDE SYDYWCEASS PHIHPVGWCK EHRRTLITPP GYPNVKHFSW DKYLEETNSL PPARAFKVK PPHGFQKKMK LEVVDKRNPM FIRVATVADT DDHRVKVHFD GWNNCYDYWI DADSPDIHPV GWCSKTGHPL QPPLSPLELM EASEHGGCST PGCKGIGHFK RARHLGPHSA ANCPYSEINL NKDRIFPDRL SGEMPPASPS FPRNKRTDAN ESSSSPEIRD QHADDVKEDF EERTESEMRT SHEARGAREE PTVQQAQRRS AVFLSFKSPI PCLPLRWEQQ SKLLPTVAGI PASKVSKWST DEVSEFIQSL PGCEEHGKVF KDEQIDGEAF LLMTQTDIVK IMSIKLGPAL KIFNSILMFK AAEKNSHNEL
Source	Baculovirus
Target Names	L3MBTL3
Protein Names	Recommended name: Lethal(3)malignant brain tumor-like protein 3 Short name= H-l(3)mbt-like protein 3 Short name= L(3)mbt-like protein 3 Alternative name(s): MBT-1
Expression Region	1-780
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.