



Recombinant Human Protein FAM84A (FAM84A)

Product Code	CSB-BP846652HU
Abbreviation	FAM84A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96KN4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGNQLDRITH LNYSELPTGD PSGIEKDEL R VGVAYFFSDD EEDLDERGQP DKFGVKAPPG CTPCPESPSR HHHHLLHQLV LNETQFSAFR GQECIFSKVS GGPQGADLSV YAVTALPALC EPGDLLELLW LQPAPEPPAP APHWAVYVGG GQIIHLHQGE IRQDSLYEAG AANVGRVNS WYRYRPLVAE LVVQNACGHL GLKSEEICWT NSESFAAWCR FGKREFKAGG EVPAGTQPPQ QQYYLKVHLG ENKVHTARFH SLEDLIREKR RIDASGRLRV LQELADLVDD KE
Source	Baculovirus
Target Names	FAM84A
Protein Names	Recommended name: Protein FAM84A Alternative name(s): Neurologic sensory protein 1 Short name= NSE1
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.