



Recombinant Human DnaJ homolog subfamily C member 30 (DNAJC30)

Product Code	CSB-MP846659HU
Abbreviation	DNAJC30
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96LL9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAMRWRWWQ RLLPWRLQA RGFPQNSAPS LGLGARTYSQ GDCSYSRTAL YDLLGVPSTA TQAQIKAAYY RQCFLYHPDR NSGSAEAAER FTRISQAYVV LGSATLRRKY DRGLLSDEDL RGPVVRPSRT PAPDPGSPRT PPPTSRTHDG SRASPGANRT MFNFDAFYQA HYGEQLERER RLRARREALR KRQEYRSMKG LRWEDTRDTA AIFLIFSIFI IIGFYI
Source	Mammalian cell
Target Names	DNAJC30
Protein Names	Recommended name: DnaJ homolog subfamily C member 30 Alternative name(s): Williams-Beuren syndrome chromosomal region 18 protein
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This intronless gene encodes a member of the DNAJ molecular chaperone homology domain-containing protein family. This gene is deleted in Williams syndrome, a multisystem developmental disorder caused by the deletion of contiguous genes at 7q11.23.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.