



Recombinant Human Neuroblastoma breakpoint family member 4 (NBPF4)

Product Code	CSB-BP846664HU
Abbreviation	NBPF4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96M43
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVVSADPLSS ERAEMNILEI NQELRSQLAE SNQQFRDLKE KFLITQATAY SLANQLKKYK CEEYKDIIDS VLRDELQSME KLAEKLRQAE ELRQYKALVH SQAKELTQLR EKLREGRDAS RVLNKHKLKTL LTPDDPKDSQ GQDLREQLAE GHRLAEHLVH KLSPENDEDE DEDEDDKDEE VEKVQESPAP REVQKTEEKE VPQDSLEECA VTCSNSHNPS NSNQPHRSTK ITFKEHEVDS ALVVESEHPH DEEEEALNIP PENQNDHEEE EGKAPVPPRH HDKSNSYRHR EVSFLALDEQ KVCSAQDVAR DYSNPKWDET SLGFLEKQSD LEEVKGQETV APRLSRGPLR VDKHEIPQES LDGCCLT PSI LPDLTPSYHP YWSTLYS FED KQVSLALVDK IKKDQEEIED QSPPCPRLS ELPEVKEQEV PEDSVNEVYL TPSVHHDVSD CHQPYSSTLS SLEDQLACSA LDVASPTEAA CPQGTWSGDL SHHQSEVQVS QAQLEPSTLV PSCLRLQLDQ GFHCGNGLAQ RGLSSTTCSF SANADSGNQW PFQELVLEPS LGMKNPPQLE DDALEGSASN TQGRQVTGRI RASLVLILKT IRRRLPFSKW RLAFRFAGPH AESAEIPNTA GRTQRMAG
Source	Baculovirus
Target Names	NBPF4
Protein Names	Recommended name: Neuroblastoma breakpoint family member 4
Expression Region	1-638
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.