



# Recombinant Human Neuroblastoma breakpoint family member 4 (NBPF4)

<b>Product Code</b>	CSB-MP846664HU
<b>Abbreviation</b>	NBPF4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96M43
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVVSADPLSS ERAEMNILEI NQELRSQLAE SNQQFRDLKE KFLITQATAY SLANQLKKYK CEEYKDIIDS VLRDELQSME KLAEKLRQAE ELRQYKALVH SQAKELTQLR EKLREGRDAS RVLNKHKLKTL LTPDDPKSQ GQDLREQLAE GHRLAEHLVH KLSPENDEDE DEDEDDKDEE VEKVQESPAP REVQKTEEKE VPQDSLEECA VTCSNSHNPS NSNQPHRSTK ITFKEHEVDS ALVVESEHPH DEEEEALNIP PENQNDHEEE EGKAPVPPRH HDKSNSYRHR EVSFLALDEQ KVCSAQDVAR DYSNPKWDET SLGFLEKQSD LEEVKGQETV APRLSRGPLR VDKHEIPQES LDGCCLTPSI LPDLTPSYHP YWSTLYSFED KQVSLALVDK IKKDQEEIED QSPPCPRLSQ ELPEVKEQEV PEDSVNEVYL TPSVHHDVSD CHQPYSSTLS SLEDQLACSA LDVASPTEAA CPQGTWSGDL SHHQSEVQVS QAQLEPSTLV PSCLRLQLDQ GFHCGNGLAQ RGLSSTTCSF SANADSGNQW PFQELVLEPS LGMKNPPQLE DDALEGSASN TQGRQVTGRI RASLVLILKT IRRRLPFSKW RLAFRFAGPH AESAEIPNTA GRTQRMAG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	NBPF4
<b>Protein Names</b>	Recommended name: Neuroblastoma breakpoint family member 4
<b>Expression Region</b>	1-638
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.