



Recombinant Human Protein HEXIM2 (HEXIM2)

Product Code	CSB-MP846669HU
Abbreviation	HEXIM2
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q96MH2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMATPNQTAC NAESPVALEE AKTSGAPGSP QTPPERHDSG GSLPLTPRME SHSEDEDLAG AVGGLGWNSR SPRTQSPGGC SAEAVLARKK HRRRPSKRKR HWRPYLELSW AEKQQRDERQ SQRASRVREE MFAKGQPVAP YNTTQFLMND RDPEEPNLDV PHGISHPGSS GESEAGDSDG RGRAHGEFQR KDFSETYERF HTESLQGRSK QELVRDYLEL EKRLSQAE TRRLQQLQAC TGQQSCRQVE ELAAEVQRLR TENQRLRQEN QMWNREGCRC DEEPGT
Source	Mammalian cell
Target Names	HEXIM2
Protein Names	Recommended name: Protein HEXIM2 Alternative name(s): Hexamethylene bis-acetamide-inducible protein 2
Expression Region	1-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>