



Recombinant Human Calcium/calmodulin-dependent protein kinase type 1G (CAMK1G)

Product Code	CSB-BP846679HU
Abbreviation	CAMK1G
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96NX5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGRKEEDDCS SWKKQTTNIR KTFIFMEVLG SGAFSEVFLV KQRLTGKLFA LKCIKKSPAF RDSSLENEIA VLKKIKHENI VTLEDIYEST THYYLVMQLV SGGELFDRIL ERGVYTEKDA SLVIQQVLSA VKYLHENGIV HRDLKPENLL YLTPEENSKI MITDFGLSKM EQNGIMSTAC GTPGYVAPEV LAQKPYSKAV DCWSIGVITY ILLCGYPPFY EETESKLFEK IKEGYEFES PFWDDISESA KDFICHLLEK DPNERYTCEK ALSHPWIDGN TALHRDIYPS VSLQIQKNFA KSKWRQAFNA AAVVHHMRKL HMNLHSPGVR PEVENRPPET QASSETSRPSS PEITITEAPV LDHSVALPAL TQLPCQHGRR PTAPGGRSLN CLVNGSLHIS SSLVPMHQGS LAAGPCGCCS SCLNIGSKGK SSYCSEPTLL KKANKKQNFK SEVMVPVKAS GSSHCragQT GVCLIM
Source	Baculovirus
Target Names	CAMK1G
Protein Names	Recommended name: Calcium/calmodulin-dependent protein kinase type 1G EC= 2.7.11.17 Alternative name(s): CaM kinase I gamma Short name= CaM kinase IG Short name= CaM-KI gamma Short name= CaMKI gamma Short name= C
Expression Region	1-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein similar to calcium/calmodulin dependent protein kinase, however, its exact function is not known.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.