



Recombinant Human Myotubularin-related protein 9 (MTMR9)

Product Code	CSB-MP847211HU
Abbreviation	MTMR9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96QG7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEFAELIKTP RVDNVVLHRP FYPAVEGTLC LTGHHLLSS RQDNTEELWL LHSNIDAIDK RFVGS LGTII IKCKDFRIIQ LDIPGMEECL NIASSIEALS TLDSITLMYP FFYRPMFEVI EDGWSHFLPE QEFELYSSAT SEWRLSYVVK EFAVCPSPYP IVTVPKSIDD EALRKVATFR HGGFRPVLSY YHKKNGMVIM RSGQPLTGTN GRRCKEDEKL INATLRAGKR GYIIDTRSLN VAQQTRAKGG GFEQEAHYPQ WRRIHKSIEK YHILQESLIK LVEACNDQTH NMDRWLSKLE ASNWLTHIKE ILTTACLAAQ CIDREGASIL IHGTEGTDST LQVTS LAQII LEPRSRITRG FEALIEREWL QAGHPFQQRC AQSAYCNTKQ KWEAPVFLLF LDCVWQILRQ FPCSFENEN FLIMLFEHAY ASQFGTFLGN NESERCKLKL QQKTMSLWSW VNQPSELSKF TNPLFEANNL VIWPSVAPQS LPLWEGIFLR WNRSSKYLDE AYEEMVNIIE YNKELQAKVN ILRRQLAELE TEDGMQESP
Source	Mammalian cell
Target Names	MTMR9
Protein Names	Recommended name: Myotubularin-related protein 9
Expression Region	1-549
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a myotubularin-related protein that is atypical to most other members of the myotubularin-related protein family because it has no dual-specificity phosphatase domain. The encoded protein contains a double-helical motif similar to the SET interaction domain, which is thought to have a role in the control of cell proliferation. In mouse, a protein similar to the encoded protein binds with MTMR7, and together they dephosphorylate phosphatidylinositol 3-phosphate and inositol 1,3-bisphosphate.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.