



Recombinant Human 1-aminocyclopropane-1-carboxylate synthase-like protein 1 (ACCS)

Product Code	CSB-EP847213HU-B
Abbreviation	ACCS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96QU6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFTLPQKDFR APTTCLGPTC MQDLGSSHGE DLEGECSRKL DQKLPELRGV GDPAMISSDT SYLSSRGRMI KFWWDSAEEG YRTYHMDEYD EDKNPSGIIN LGTSSENKLCF DLLSWRLSQR DMQRVEPSLL QYADWRGHLF LREEVAKFLS FYCKSPVPLR PENVVVLNGG ASLFSALATV LCEAGEAFLI PTPYYGAIQ HVCLYGNIRL AYYVLDSEVT GLDTRPFQLT VEKLEMALRE AHSEGKVKVG LILISPNPL GDVYSPEELQ EYLVFAKRHR LHVIVDEVYM LSVFEKSVGY RSVLSLERLP DPQRTHVMWA TSKDFGMSGL RFGTLYTENQ DVATAVASLC RYHGLSGLVQ YQMAQLLRDR DWINQVYLPE NHARLKAHT YVSEELRALG IPFLSRGAGF FIWVDLRKYL PKGTFFEEML LWRRFLDNKV LLSFGKAFEC KEPGWFRFVF SDQVHRLCLG MQRVQQVLAG KSQVAEDPRP SQQEPSDQR R
Source	E.coli
Target Names	ACCS
Protein Names	Recommended name: 1-aminocyclopropane-1-carboxylate synthase-like protein 1 Short name= ACC synthase-like protein 1
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.