



Recombinant Human Cyclin-L2 (CCNL2)

Product Code	CSB-EP847223HU-B
Abbreviation	CCNL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96S94
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p> AAAAAAAGA AGSAAPAAAA GAPGSGGAPS GSQGV LIGDR LYSGVLITL E NCLLPDDKLR FTPSMSSGLD TDTETDLRVV GCELIQAAGI LLRLPQVAMA TGQVLFQRFF YTKSFVKHSM EHVSMACVHL ASKIEEAPRR IRDVINVFHR LRQLRDKKKP VPLLLDQDYV NLKNQIIKAE RRVLKELGFC VHVKHPHKII VMYLQVLECE RNQHLVQTSW NYMNDLRTD VVRFQPE SI ACACIYLAAR TLEIPLPNRP HWFLLFGATE EEIQEICLKI LQLYARKKVD LTHLEGEVEK RKHAIEEAKA QARGLLPGGT QVLDGTSGFS PAKLVESPK EGKGSKPSPL SVKNTKRRLE GAKKAKADSP VNGLPKGRES RSRRSRSREQS YSRSPRSAS PKRRKSDSGS TSGGSKSQSR SRSRSDSPPR QAPRSAPYKG SEIRGSRKSK DCKYPQKPHK SRSRSSSR SR SRERADNP GK YKKKSHYY RDQRRERSRS YERTGRRYER DHPGHSRHR </p>
Source	E.coli
Target Names	CCNL2
Protein Names	Recommended name: Cyclin-L2 Alternative name(s): Paneth cell-enhanced expression protein
Expression Region	2-520
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.