



# Recombinant Human Spermatid perinuclear RNA-binding protein (STRBP)

<b>Product Code</b>	CSB-BP847225HU
<b>Abbreviation</b>	STRBP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96SI9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MRSIRSFAND DRHVMVKHST IYPSPEELEA VQNMVSTVEC ALKHVSDWLD ETNKGTKTEG ETEVKKDEAG ENYSKDQGGG TLCGVMRIGL VAKGLLIKDD MDLELVLMCK DKPTETLLNT VKDNLPIQIQ KLTEEKYQVE QCVNEASIII RNTKEPTLTL KVILTSPLIR DELEKKDGEN VSMKDPPDLL DRQKCLNALA SLRHAKWFQA RANGLKSCVI VLRILRDLN RVPTWAPLKG WPLELICEKS IGTCNRPLGA GEALRRVMC LASGILLPGG PGLHDPCERD PTDALSYMTI QQKEDITHSA QHALRLSAFG QIYKVLEMDP LPSSKPFQKY SWSVTDKEGA GSSALKRPF E DGLGDDKDPN KKMKNRNLKI LDSKAIDL MN ALMRLNQIRP GLQYKLLSQS GPVHAPVFTM SVDVDGTTYE ASGPSKKTAK LHVAVKVLQA MGYPTGFDAD IECMSSDEKS DNESKNETVS SNSSNNTGNS TTETSSTLEV RTQGPILTAS GKNPVMELNE KRRGLKYELI SETGGSHDKR FVMEVEVDGQ KFRGAGPNKK VAKASAALAA LEKLFSGPNA ANNKKKKIIP QAKGVVNTAV SAAVQAVRGR GRGTLTRGAF VGATAAPGYI APGYGTPYGY STAAPAYGLP KRMVLLPVMK FPTYPVPHYS FF
<b>Source</b>	Baculovirus
<b>Target Names</b>	STRBP
<b>Protein Names</b>	Recommended name: Spermatid perinuclear RNA-binding protein
<b>Expression Region</b>	1-672
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .