



Recombinant Human UPF0464 protein C15orf44 (C15orf44)

Product Code	CSB-MP847227HU
Abbreviation	C15orf44
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96SY0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPTVVVMDVS LSMTRPVSIE GSEEYQRKHL AAHGLTMLFE HMATNYKLEF TALVVFSSLW ELMVPFTRDY NTLQEALSNM DDYDKTCLES ALVGVCNIVQ QEWGGAIPCQ VVLVTDGCLG IGRGSLRHSL ATQNQRSESN RFPLPPFPFS KLYIMCMANL EELQSTDSLE CLERLIDLNN GEGQIFTIDG PLCLKNVQSM FGKLIDLAYT PFHAVLKCGH LTADVQVFPR PEPFVVDEEI DPIPKVINTD LEIVGFIDIA DISSPPVLSR HLVLPIALNK EGDEVGTGIT DDNEDENSAN QIAGKIPNFC VLLHGSLKVE GMVAIVQLGP EWHGMLYSQA DSKKKS NLMM SLFEPGPEPL PWLGKMAQLG PISDAKENPY GEDDNKSPFP LQPKNKRSYA QNVTVWIKPS GLQTDVQKIL RNARKLPEKT QTFYKELNRL RKAALAFGFL DLLKGVADML ERECTLLPET AHPDAAFQLT HAAQQLKLAS TGTSEYAAAYD QNITPLHTDF SGSSTERI
Source	Mammalian cell
Target Names	INTS14
Protein Names	Recommended name: UPF0464 protein C15orf44
Expression Region	1-518
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.